



Sealed Air

Cryovac, Inc.
Cryovac Food Packaging Division
100 Rogers Bridge Road
P.O. Box 464
Duncan, SC 29334-0464

Tel: 864-433-2000
Fax: 864-433-3146

October 27, 2005

Mr. Dennis O'Donnell
FRONTIER PACKAGING
1201 Andover Park East
Seattle, WA 98188

Dear Mr. O'Donnell:

Please be advised that the Cryovac HP2700 (BDX 9124) bags shipped to your firm have a measured Oxygen Transmission Rate (OTR) greater than 10,000 cc (STP)/(24 hours, m², atm) according to guidelines listed in ASTM D3985 (TSR#28869).

This OTR complies with the FDA guidance published in *Fish and Fisheries Products Hazard and Control Guidance* (Third Edition, June 2001), and the compliance policies stated in FDA Import Alert #16-125 (9/25/02).

All bags will be printed with the "CRYOVAC 10K OTR" logo on the backside of the bag. This will alert state and federal inspectors at any point in the retail, food service, or import distribution channels that the package permeability meets FDA guidelines. Documentation from Cryovac® should be maintained at your facility that certifies that all packaging materials with the "CRYOVAC 10K OTR" logo have a measured OTR greater than 10,000 cc/m²/24 hr @STP.

Please let me know if you have any questions or need additional information.

Sincerely yours,

Clifford A. Kendrick
Seafood Marketing & Applications Specialist